


# AnyOne Surgical Kit

## I. AnyOne Internal / External / OneStage Kit

Ref.C

KA0IN3003



## PDF Compressor Free Version

### Shaping Drill


- Each drill has depth marking lines from 7.0mm to 15.0mm
- The dual marking system (grooves and laser markings) provides visual and radio graphic depth verification during surgery.



Drill Diameter	Ø2.8	Ø3.3	Ø3.6	Ø4.2	Ø4.8	Ø5.8	Ø6.9
Y length	0.58	0.59	0.68	0.85	0.89	0.94	0.94

※ Actual drill length : Drill length does not normally include the Y dimension of the drill.

※ Markings on the Shapping Drill are 0.5mm longer than the fixture so fixtures will automatically be placed 0.5mm subcrestally if the drilling protocol is followed.




※ To place a Ø5.0 x 10mm length fixture, the required bone depth would be 10.89mm.

For example : 0.5mm(subcrestal concept) + 0.89mm(Y dimension of drill tip) + 9.5mm (fixture length)

### Dense Drill


- To control initial stability in dense bone (type I & II), use the Dense Drill to remove and shape the cortical bone.



## PDF Compressor Free Version

### ► ► Surgical drilling sequence


- AnyOne fixtures offer optimum initial stability when they are used with the following drill sequence guide, AnyOne implants should be placed 0.5mm subcrestally.



#### Ø3.5 Fixture

##### Ø3.5 drilling sequence


10.0mm is the fixture length, The Shaping Drills are 0.59mm longer than the fixture, so total drill depth is 10.59mm.



#### Ø4.0 Fixture

##### Ø4.0 drilling sequence

10.0mm is the fixture length, The Shaping Drills are 0.68mm longer than the fixture, so total drill depth is 10.68mm.




#### Ø4.5 Fixture

##### Ø4.5 drilling sequence

10.0mm is the fixture length, The Shaping Drills are 0.85mm longer than the fixture, so total drill depth is 10.85mm.


## PDF Compressor Free Version



### Ø5.0 Fixture

#### Ø5.0 drilling sequence


10.0mm is the fixture length, The Shaping Drills are 0.89mm longer than the fixture, so total drill depth is 10.89mm.



### Ø6.0 Fixture

#### Ø6.0 drilling sequence

10.0mm is the fixture length, The Shaping Drills are 0.94mm longer than the fixture, so total drill depth is 10.94mm.



### Ø7.0 Fixture



#### Ø7.0 drilling sequence

10.0mm is the fixture length, The Shaping Drills are 0.94mm longer than the fixture, so total drill depth is 10.94mm.

## PDF Compressor Free Version


### ►► Handpiece & Ratchet Connector

#### Internal





- \* Fixture level : Placement should be 0.5mm subcrestal.
- \* Gingival level : The line is 3mm above the bone level line and 3.5mm above the platform line

#### External



- \* Fixture Level : Place in accordance with bone level
- \* Gingival Level : Place 3mm above the bone level

#### OneStage



- \* Fixture Level : Place in accordance with bone level
- \* Gingival Level : Place 3mm above the bone level

## PDF Compressor Free Version

### ➔ Surgical Kit Components (Continued)

#### Initial Drill

- Used to pierce the cortical bone initially.
- Advisable to go into the bone to the full length of a fixture.

Diameter	Length(mm)	Ref.C
Ø1.8	33	ID1818S
	38	*ID1818M
	43	*ID1818L

(\*) Separate sales item.



#### Shaping Drill

- Each drill has depth marking lines from 7.0mm to 15.0mm.
- The dual marking system(grooves and laser markings) provides visual and radiographic depth verification during surgery.
- TiN coating on drills : Enhanced corrosion resistance and abrasion resistance.

Diameter	Length(mm)	Ref.C
Ø2.0	33	SD2018S
	38	*SD2018M
	43	*SD2018L
Ø2.5	33	SD2518S
	38	*SD2518M
	43	*SD2518L
Ø2.8	33	SD2818S
	38	*SD2818M
	43	*SD2818L
Ø3.3	33	SD3318S
	38	*SD3318M
	43	*SD3318L
Ø3.6	33	SD3618S
	38	*SD3618M
	43	*SD3618L
Ø4.2	33	SD4218S
	38	*SD4218M
	43	*SD4218L
Ø4.8	33	SD4818S
	38	*SD4818M
	43	*SD4818L
Ø5.8	33	SD5818S
	38	*SD5818M
	43	*SD5818L
Ø6.9	33	SD6918S
	38	*SD6918M
	43	*SD6918L

(\*) Separate sales item.



#### Dense Drill

- Used to remove and shape cortical bone to control initial stability in dense bone (type I & II).
- TiN coating on drills : Enhanced corrosion resistance and abrasion resistance.

Diameter	Type	Ref.C
Ø3.9	Long	DD39
Ø4.3		DD43
Ø4.8		DD48
Ø5.3		DD53
Ø6.3		DD63
Ø7.3		DD73



## PDF Compressor Free Version

### ➔ Surgical Kit Components (Continued)

#### Handpiece Connector

- Used with Handpiece to remove fixture from ampule and to place the fixture.
- Spring type connection allows for easy and secure pick-up and positioning of the fixture.
- First mark on the shaft indicate the position of the fixture platform, For second mark, the bottom of the black line is 3mm and the top of the black line is 4mm(from fixture platform).
- Especially useful in flapless surgery.

##### AnyOne Internal& External

Length(mm)	Type	Connection	Ref.C
5	*Ultra-short	Hex. 2.5	HCU25
10	Short		HCS25
15	Long		HCL25

(\*) Separate sales item




##### MiNi

Length(mm)	Type	Connection	Ref.C
10	Short	Hex. 1.7	HCS17
15	Long		HCL17

##### OneStage

Length(mm)	Type	Connection	Ref.C
6	*Ultra-short	Octa. 3.1	MTHC200U
9	Short		MTHC200S
16	Long		MTHC200L

(\*) Separate sales item



#### Ratchet Connector

- Used for inserting or removing a fixture with the Ratchet Wrench.
- Check to make sure the Ratchet Connector is completely seated into the Ratchet Wrench before using.
- Excessive force can cause damage to internal hex of fixture.
- Marks on the shaft indicate the position of fixture platform. Bottom of the black line is 3mm and top of black line is 4mm(from fixture platform).
- Especially useful in flapless surgery.


##### Internal& External

Length(mm)	Type	Connection	Ref.C
10	*Ultra-short	Hex. 2.5	RCU25
15	Short		RCS25
20	Long		RCL25

(\*) Separate sales item

##### MiNi

Length(mm)	Type	Connection	Ref.C
15	Short	Hex. 1.7	RCS17
20	Long		RCL17




#### Final Driver

- Used to attach or remove the fixture by connecting to Ratchet Wrench
- Used to mount the Ratchet Connector fully on the Ratchet Wrench

##### OneStage

Length(mm)	Type	Connection	Ref.C
6	Short	Octa. 3.1	MOHD310S
13	Long		MOHD310





PDF Compressor Free Version

Hand Driver (1.2 Hex)

- Used for all Cover Screws, Abutment Screws, and Healing Abutments.
- Available in 4 lengths for added convenience.
- Hand Driver can be directly inserted into the Torque Wrench without using an adaptor.
- Hex tip can with stand 35-45Ncm of torque without distorting.

Length(mm)	Type	Ref.C
5	*Ultra-short	TCMHDU1200
10	Short	TCMHDS1200
15	Long	TCMHDL1200
20	*Extra-long	TCMHDE1200

(\*) Separate sales item



Hand Driver (0.9 Hex)

- Used for AnyOne External fixture cover screw.
- Available in 3 lengths for convenience.
- Hand Driver can be directly inserted in the to Torque Wrench without using an adaptor.
- Hex tip can with stand 25-35Ncm of torque without distorting.

Length(mm)	Type	Ref.C
5	*Ultra-short	TCMHDU0900
10	Short	TCMHDS0900
15	Long	TCMHDL0900

(\*) Separate sales item



Drill Extension

- No more than 35Ncm torque : May distorted when excessive force is applied.
- Extends drills & other handpiece instruments.

Ref.C
MDE150






## PDF Compressor Free Version

### ➔ Surgical Kit Components

#### Direction Indicator

- Confirms drilling direction and functions as a parallel guide for additional osteotomies.
- Each end of the Direction Indicator has a different diameter
  - Ø2.0 and Ø2.8.


Diameter	Ref.C
Ø2.0 / Ø2.8	MDI100



#### Path Finder

- After the fixture is placed, a Path Finder may be connected into the fixture and function as a parallel guide for additional osteotomies.
- Grooves indicate the distance from the fixture platform. The first groove is 0.3mm and the second groove is 1mm, especially useful in flapless surgery.

Length(mm)	Ref.C
15	PF



#### Ratchet Wrench

- Used to exert more force than the Handpiece.
- No bearing system : No breakage and no corrosion problems.
- Arrow laser marking indicates direction of force.

Ref.C
MRW040S

