

# **R2GATE Surgical KIT**







#### 96 R2GATE FULL Surgical KI

- 496 I. AnyRidge System
- 498 II. AnyOne System

#### 506 R2GATE Standard KIT

- 506 I. TSIII System (Osstem co.)
- 508 II. SuperLine System (Dentium co.)
- 510 III. UFII System (DIO co.)
- 512 IV. ISII System (Neo Biotech co.)

#### 520 R2GATE Universal KIT

- 524 1. AnyRidge Octa 1
- 524 **2. Straumann**
- 526 3. Nobel Biocare
- 526 **4. Astra**
- 527 **5. Biomet 3i**
- 527 **6. TSIII**
- 28 7. SuperLine
- 528 **8. ISII**
- <sup>529</sup> 9. UFII
- 530 **10. Final Drill Option**


#### 532 R2GATE Narrow KIT


34 R2GATE Anchor KIT


## R2GATEmpFullrSurgical KIT

# I. R2GATE Full Surgical Kit for AnyRidge System

- If you only use a specific system, corresponding system's full kit can be provided.
- R2GATE full surgical kit is composed with all of drills and components that are needed for the Digital Guided Surgery which uses R2GATE Guide™ after R2GATE™ diagnosis. It helps to actualize minimally invasive surgery and makes exact clinical result as the diagnosis.











## **II. R2GATE Full Surgical Kit** for AnyOne System

- · If you only use a specific system, corresponding system's full kit can be provided.
- nents that are needed for the Digital Guided Surgery which uses minimally invasive surgery and makes exact clinical result as the



**Cortical Bone Drill** 




## Components for R2GATE Full Surgical Kit (Continued)

- If you only use a specific system, corresponding system's full kit can be provided.
- R2GATE full surgical kit is composed with all of drills and components that are needed for the Digital Guided Surgery which uses R2GATE Guide<sup>TM</sup> after R2GATE<sup>TM</sup> diagnosis. It helps to actualize minimally invasive surgery and makes exact clinical result as the diagnosis.

#### **Initial Drill**

- Use the initial drill in order to mark the drilling position on the bone. Start drilling slowly, when drill guide part is fully contacted with drilling core of R2GATE Guide™.
- Recommended drilling speed range is 300 ~ 800 RPM with copious irrigation.

Diameter	Guide Diameter	Length(mm)	Ref.C
Ø2.6	Ø5.0	1.0	R2ID2601



#### Second Drill

- This unique step-drill(from  $\emptyset$ 2.0 to  $\emptyset$ 4.6) is used to flare out the upper cortical bone of the osseotomy.
- It helps not only the rest drilling procedure but abut- ment connection. In case of hard bone, if the 2nd drilling will be disturbed by thick cortical bone. Stop the drilling and try it after final drilling procedure.

Diameter	Guide Diameter	Length(mm)	Ref.C
Ø2.5	Ø5.0	5.0	R2SD2505



### Stopper Drill

- Universal drills consist of Ø2.0, Ø.2.5, Ø2.8 diameter to enlarge the osteotomy gradually.
- The length of drill are designed as 7.0, 8.5, 10, 11.5,13mm for most common length of implant system.
- Recommended drilling speed range is 500 ~ 800 RPM with copious irrigation.

Diameter	Guide Diameter	Length(mm)	Ref.C
		6.5	AGSD2007
		8.0	AGSD2008
Ø2.0		9.5	AGSD2010
		11.0	AGSD2011
		12.5	AGSD2013
		6.5	AGSD2507
	Ø5.0	8.0	AGSD2508
Ø2.5		9.5	AGSD2510
		11.0	AGSD2511
		12.5	AGSD2513
		6.5	AGSD2807
		8.0	AGSD2808
Ø2.8		9.5	AGSD2810
		11.0	AGSD2811
		12.5	AGSD2813



#### **Bone Profiler**


• Recommended drilling speed is 300 ~ 800 RPM.

Diameter	Guide Diameter	Ref.C	
Ø4.0	Ø5.0	AGBP40	
Ø5.0		AGBP50	
Ø6.0	Ø6.5	AGBP60	



## Stopper Drill[AR]

• Recommended drilling speed is 300 ~ 800 RPM.



Diameter	Guide Diameter	Length(mm)	Ref.C
		6.5	ARSD3307
		8.0	ARSD3308
Ø3.3		9.5	ARSD3310
		11.0	ARSD3311
	Ø3.8 Ø5.0	12.5	ARSD3313
		6.5	ARSD3807
		8.0	ARSD3808
Ø3.8		9.5	ARSD3810
		11.0	ARSD3811
		12.5	ARSD3813
		6.5	ARSD4307
		8.0	ARSD4308
Ø4.3		9.5	ARSD4310
		11.0	ARSD4311
		12.5	ARSD4313

Diameter	Guide Diameter	Length(mm)	Ref.C
		6.5	ARSD4807
		8.0	ARSD4808
Ø4.8		9.5	ARSD4810
		11.0	ARSD4811
		12.5	ARSD4813
	Ø6.5	6.5	ARSD5407
		8.0	ARSD5408
Ø5.4		9.5	ARSD5410
		11.0	ARSD5411
		12.5	ARSD5413
		6.5	ARSD5908
		8.0	ARSD5907
Ø5.9		9.5	ARSD5910
		11.0	ARSD5911
		12.5	ARSD5913

### Stopper Drill[AO]

• Recommended drilling speed is 300  $\sim$  800 RPM.



Diameter	Guide Diameter	Length(mm)	Ref.C
		7.0	AOSD3307
		8.0	AOSD3308
Ø3.3		9.5.0	AOSD3310
		11.0	AOSD3311
		12.5	AOSD3313
		7.0	AOSD3607
	Ø5.0	8.0	AOSD3608
Ø3.6		9.5	AOSD3610
		11.0	AOSD3611
		12.5	AOSD3613
		7.0	AOSD4207
		8.0	AOSD4208
Ø4.2		9.5	AOSD4210
		11.0	AOSD4211
		12.5	AOSD4213


Diameter	Diameter	Length(mm)	Ref.C
		7.0	AOSD4807
		8.0	AOSD4808
Ø4.8		9.5	AOSD4810
		11.0	AOSD4811
	005	12.5	AOSD4813
	Ø6.5	7.0	AOSD5807
		8.0	AOSD5808
Ø5.8		9.5	AOSD5810
		11.0	AOSD5811
	12.5	AOSD5813	

## Components for R2GATE Full Surgical Kit (Continued)

#### Cortical Bone Drill[AR]

• Recommended drilling speed : 300 ~ 800 RPM

Diameter	Guide Diameter	Length(mm)	Ref.C
Ø3.4			R2CD3405
Ø3.8	Ø5.0		R2CD3805
Ø4.3	25.0	5.0	R2CD4305
Ø4.8			R2CD4805
Ø5.3			R2CD5305
Ø5.8	Ø6.5		R2CD5805
Ø6.3			R2CD6305



#### Cortical Bone Drill[AO]

• Recommended drilling speed : 300 ~ 800 RPM

Diameter	Guide Diameter	Length(mm)	Ref.C
Ø3.9			AODD39
Ø4.3	Ø5.0	6.0	AODD43
Ø4.8			AODD48
Ø5.3	Ø6.5	5.5	AODD53
Ø6.3			AODD63



#### Implant Carrier[AR]

- The purpose of tab drills in the universal kit system is insertion test, some of implant are required this procedure before final fixture insertion, choose the one-step under size of tab to protect from enlargement of osteotomy.
- Recommended insertion torque and speed is 45  $\sim$  50Ncm, under 40 RPM.

Connection	Guide Diameter	Туре	Ref.C
2.3 Hex	Ø5.0	Ratchet	ICRH2324
	Ø6.5		ICWH2324
	Ø5.0	Handpiece	ICRH2324H
	Ø6.5		ICWH2324H



#### Implant Carrier[AO]

- Two different implant carriers for regular stent since Ø3.5 fixture has different abut- ment connection
- To pick up the fixture from the ampule and insert it to the ossetomy. Then turn it to clock-wise direction 2~3 times manualy.
- When it gets fixation from the osteotomy, connect the handpiece adaptor and use implant motor.
- Recommended insertion torque is 45~50Ncm.

Connection	Guide Diameter	Туре	Ref.C
	Ø5.0	Ratchet	ICRH2518
			ICRH2523
2.3 Hex	Ø6.5		ICWH2523
2.3 NeX	Ø5.0	Handpiece	ICRH2518H
			ICRH2523H
	Ø6.5		ICWH2523H



#### Carrier-Handpiece Adapter

• Useful to use the handpiece for the implant placement following initial delivery of a fixture with a fixture carrier.

Diameter	Ref.C	
5.0	AGHA	



#### **Carrier Extension**

• To extend the length of implant carrier.

Diameter	Ref.C
4.0	MRE400S



#### **Drill Extension**

- No more than 35Ncm torque : May distorted when excessive force is applied.
- Extends drills & other handpiece instruments.






#### Hand Driver (1.2 Hex)

- Used for all Cover Screws, Abutment Screws, and Healing Abutments.
- Available in 4 lengths for added convenience.
- Hand Driver can be directly inserted into the Torque Wrench without using an adaptor.
- Hex tip can with stand 35-45Ncm of torque without distorting.

Length(mm)	Туре	Ref.C
5.0	*Ultra-short	TCMHDU1200
10	Short	TCMHDS1200
15	Long	TCMHDL1200
20	*Extra-long	TCMHDE1200





#### Ratchet Wrench

- Used to exert more force than the Handpiece.
- No bearing system : No breakage and no corrosion problems.
- · Arrow laser marking indicates direction of force.





## System Options for the AnyOne External

#### Cortical Bone Drill

• Recommended drilling speed: 300~800 rpm

Diameter	Guide Diameter	Length(mm)	Ref.C
Ø3.9		Ø5.0 6.0	AODD39
Ø4.3	Ø5.0		AODD43
Ø4.8			AODD48
Ø5.3	Ø6.5	5.5	AODD53
Ø6.3			AODD63



#### Tap Drill

- The purpose of tab drills in the R2GATE universal kit system is insertion test. some of implant are required this procedure before final fixture insertion. choose the one-step under size of tab to protect from enlargement of osteotomy.
- Recommended insertion torque and speed is 45 ~50Ncm, under 40 RPM.

Diameter	Guide Diameter	Length(mm)	Ref.C
Ø3.9		9.5	R2TD35AO
Ø4.3	Ø5.0		R2TD40AO
Ø4.8			R2TD45AO
Ø5.3	Ø6.5		R2TD50AO
Ø6.3	20.5		R2TD60AO



#### Implant Carrier

- To pick up the fixture from the ampule and insert it to the ossetomy. Then turn it to clock-wise direction 2~3 times manualy.
- When it gets fixation from the osteotomy, connect the handpiece adaptor and use implant motor
- Recommended insertion torque is 45~50Ncm.

Diameter	Guide Diameter	Туре	Ref.C
0.E.Hov	Ø5.0	Ratchet	ICRAOE
2.5 Hex	Ø6.5		ICWAOE

