


# Abutment & Prosthetic Options

## I. Fixture Level Prosthesis



## PDF Compressor Free Version

### ➔ Abutment Options (Continued)

#### Lab Analog

- Replicates the fixture.
- Magenta analog for Ø3.5 fixture.
- Blue analog for all fixture sizes for Ø4.0~Ø8.0.

Type	Color	Ref.C
Small	Magenta	LA350H
Regular & Wide	Blue	LA400H




#### Temporary Abutment (Titanium)

- Abutment Screw(AS20) included

- For making provisional restoration.
- Available in both hex and non-hex.
- Grooved surface on abutment post allows for better retention of resin or wax.
- Recommend torque : 25Ncm

Profile Diameter	Height(mm)	Type	Ref.C
Ø4.5	11.0	Hex	TA4511HT
		Non-Hex	TA4511NT




#### Temporary Abutment (POM)

- Abutment Screw(AS20) included

- For making chairside provisionals for the aesthetic zone.
- Especially useful for immediate placement after extraction.
- Available in both hex and non-hex.
- Recommend torque : 25Ncm

Profile Diameter	Height(mm)	Type	Ref.C
Ø4.5	11.0	Hex	TA4511HPT
		Non-Hex	TA4511NPT



#### Fuse Abutment

- Abutment Screw(AS20)+Fuse cap included.

- For the design concept and rationale of the Fuse Abutment, please refer to Page 078.
- Recommend torque : 25Ncm

Diameter	C:H (mm)	P:H (mm)	Type	Ref.C
Ø5.5	Ø5.5	5.5	Straight	AOFAP5545P
	Ø4.5	7	15°	AOFAA5415P
	Ø4.5	7	25°	AOFAP5425P

NEW : 4mm cuff height available  
→ Adequate for deeply placed implants or thick gingival cases

